

## Hill-Rom®

### Art of Care® 4-Position Recliner



1 As shown

*R4SP: Fabric; Silverado Ginger (L74)  
Finish Honey Maple*



2



3



4

#### At a glance

- Four Positions – seating (1), TV (2), recline (3) and feet above head (4)
- Exceptionally smooth and easy to recline mechanism
- Front locking casters, positioning handles and arm grips to facilitate egress/ingress
- Accessory hooks on both sides
- SteelLock™ construction ensures superior structural integrity
- Optional removable arm, tray table or accessory bracket on left, right or both
- Choose from either wood or urethane arms

#### Hill-Rom's 4-Position Recliner delivers style, comfort and functionality for patients and caregivers.

Created to provide an aesthetically pleasing design that will enhance the overall experience for the patient and visitor. Hill-Rom's 4-Position recliner integrates with our four Art of Care® styles. Using spring box construction for seating comfort and an exceptionally smooth and easy-to-use mechanism. This recliner provides performance, durability and utility that you have come to expect from Hill-Rom. Positioning handle, front locking casters, accessory hooks, optional removable arm, tray tables, urethane arms caps and accessory brackets for an IV pole or O<sub>2</sub> tank are just a few of the distinctive options.

**Hill-Rom**

Enhancing Outcomes for Patients and Their Caregivers.™

# Hill-Rom®

## Art of Care® 4-Position Recliner

Hill-Rom's 4-Position Recliner has standard features or options that will ensure a comfortable, functional and stylish sitting surface.



### Accessory Hooks

Located on both sides, the spring loaded hooks retract automatically in place for cleanliness and safety.



### Smooth and Protected Mechanism

User friendly and smooth, the 4-position recliner mechanism is built for durability. For increased safety, the integrated ottoman mechanism is protected.



### Arm Grip

Extended arm grips assist user during ingress and egress.



### Optional Accessory Bracket

For IV pole or O<sub>2</sub> tank.



### Optional Tray Table

Conveniently placed on any side the tray table rotates easily out of the way to be stored.



### Optional Removable Arm

Located on either side or both, the removable arm allows for ingress and egress.

Product	CAL 133	Options	Caps Finish	Upholstery
<b>PLEASE SPECIFY THE FOLLOWING FOR EACH 4 POSITION RECLINER</b>				
R4SP	Cal 133 available	Right or left tray table. Right, left or both removable arm. Right, left or both accessory bracket. Urethane arms.	See Hill-Rom. finish options brochure or visit <a href="http://www.hill-rom.com">www.hill-rom.com</a>	See Hill-Rom. fabric cards or visit <a href="http://www.hill-rom.com">www.hill-rom.com</a>

Visit [www.hill-rom.com](http://www.hill-rom.com) to "Design Your Own" Art of Care® 4-Position Recliner

## Technical Specifications



### R4SP 4-Position Recliner

W 29.00 D 40.00 H 50.00 SW 23.00 SD 22.00 SH 19.00

### Footprint Fully Reclined

W 29.00 D 69.00

COM Ydg: 6.25

Maximum Weight: 300 lbs

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2009 Hill-Rom Services, Inc ALL RIGHTS RESERVED.

151602 rev 2 8/17/09



Enhancing Outcomes for Patients and Their Caregivers.™

USA 800-445-3730  
Canada 800-267-2337

[www.hill-rom.com](http://www.hill-rom.com)