



samaritan® AED

Price List - June 2005

Delivery: FOB Factory San Clemente, CA, USA.

Part No.	Description	Price U.S.D.
Samaritan® AED Defibrillators		
Samaritan® AEDs incorporate a lightweight (4.2 lbs), durable-rubberized design with a low energy, escalating SCOPE Biphasic waveform.		
SAM 001 - **	Samaritan® AED - with ECG / Text Display & Manual Override Options	\$2,695.00
SAM 002 - **	Samaritan® AED - with ECG / Text Display Option	\$2,595.00
SAM 003 - **	Samaritan® AED - with Text only Display Option	\$2,495.00
Samaritan® AEDs include user manual, 2 pairs defibrillation pads, Data-Pak Battery pack, quick reference instruction card and standard 5 Year Warranty.		
Samaritan® Data Pak Battery		
SBP 001	Samaritan® Data-Pak™ Battery Pack w/32 MB Internal Memory provides over 12 hours of continuous monitoring or 120 shocks. 5 Year Shelf life, 2 Year Stand-by life. Integral memory system provides over 12 hours of ECG and event recording and 60 minutes of audio recording. Also serves as device configuration module.	\$245.00
SBP 002	Samaritan® Data-Pak™ + Rechargeable Battery Provides at least 5 hours of continuous monitoring, or 70 shocks at 200J. 6 month stand-by life. Recharge time 2 hours or less.	\$275.00
SBP 010	Battery Charger for Samaritan® Data-Pak™	\$115.00
Samaritan® Defibrillation Pads		
SDE 201	Samaritan® Defibrillation Pads - 1 set of pads	\$25.00
SDE 210	Samaritan® Defibrillation Pads - 10 sets of pads	\$250.00
Samaritan® Data Management		
SDU 001	Samaritan® Data Recovery Unit (SDRU) & SAVER™ Software System	\$395.00
SDU 002	Additional SAVER™ Data Management & Review Software Licenses	\$395.00
SDU 003	Additional Samaritan® Data Recovery Units (SDRU)	\$145.00
Samaritan® Carrying Cases & Storage Systems		
SSS 001	Soft Carrying Case - Samaritan AED	\$95.00
SSS 002	Hard Carrying Case - Samaritan AED	\$245.00
SSS 011	Rescue Cabinet with Alarm - Deluxe (Designed to fit AED & O2 tank)	\$465.00
SSS 012	Rescue Cabinet with Alarm - Basic (Designed to fit AED)	\$295.00
SSS 013	Wall Bracket - Samaritan® AED	\$125.00
SSS 014	Vehicle Bracket - Samaritan® AED	\$155.00
SSS 021	AED Wall Sign - includes triangular wall mount	\$49.00
Samaritan® Accessories		
SAC 011	AED Prep Kit (Includes: gloves, barrier mask, scissors, razor, antimicrobial Wipe, absorbent towel)	\$29.00
SAC 010	3-Lead ECG Monitor Adaptor - IEC Code	\$155.00
Samaritan® Training System		
STU 001	Samaritan® AED Trainer - realistic Samaritan® AED appearance complete with 13 training scenarios, rechargeable battery & charger, trainer defibrillator pads, user guide and carrying case.	\$445.00
STU 002	Samaritan® Trainer replacement battery charger	\$129.00
STU 003	Samaritan® Trainer defibrillation pad (1 set)	\$20.00
STU 004	Samaritan® Trainer Nylon Carry Case	\$65.00
STU 005	Samaritan® ECG / Defibrillation Simulator - provides 5 ECG waveforms (VF, VT Fast, VT Slow, NSR, Asystole) with auto-convert feature, Low Battery and Shock indicator.	\$445.00