

Hill-Rom Basic Care™ Hospital Beds



At a glance

- Made to Hill-Rom's high quality standards for design, engineering, and manufacturing
- Designed and constructed to meet stringent international quality safety standards
- Manual and electric versions available
- Standard length and Special length
- Battery backup (electric version)
- Patent pending on innovative siderails
- Caregiver and patient controls (electric version)
- CPR release (electric version)

Hill-Rom Quality and Safety at an Affordable Price.

Hill-Rom created the Basic Care™ Bed to deliver Hill-Rom quality and safety at an affordable price, because you and your patients deserve the highest quality in patient care equipment. The Basic Care™ Bed meets stringent international quality and safety standards, and contains important safety features, to assist the caregiver in providing superior patient care.

Hill-Rom Basic Care™ Hospital Beds



Clocking, TuckAway™,
Zero-transfer Gap Siderails

Instant CPR Release
(electric version)



Battery Backup

Four-Position
IV Pole Mounts



Siderail Caregiver Controls
(electric only)

Central Brake and
Steer (optional)





Siderail Patient and Caregiver Controls



Electric Foot End Caregiver Controls – Standard interface controls make it easy for caregivers to operate the bed. (electric version)



Blow-molded Headboard and Footboard with Single-Step Removal



5" (12.7 cm) Casters



Manual Foot End Caregiver Controls



3 DC Motors (electric version)





Standard Features

- Clocking, TuckAway™, Zero-transfer gap, siderails
- Battery back-up (electric version)
- Non-powered Instant CPR Release (electric version)
- Caregiver and patient siderail controls (electric version)
- Caregiver foot end controls
- Complete bed articulation: Hi-low, head, knee, automatic contour (automatic contour for electric version)
- Caregiver articulation lockout controls (electric version)
- 4-caster brake and steer caster
- 12° Trendelenburg and reverse-Trendelenburg
- Articulation indicator
- Four-position IV pole mounts
- Blow-molded headboard and footboard, single-step removal
- Headboard and footboard push handles
- 3 DC motors (electric version)
- Mobile (5"/12.7 cm) casters
- Corner roller bumpers
- Standard length sleep deck - 80" (203 cm)
- Compatibility with most:
 - Standard Mattresses
 - Overbed tables in the low position
 - IV Poles
 - Sling scale devices
 - Hoists
 - Foley-type bags
 - Urinal devices

Optional Features

- Central brake and steer
- Special length sleep deck – 76" (193 cm)
- IV Pole

Meets these internationally recognized quality and safety standards:

- **IEC 60601-1** – Medical Electrical Equipment: General Requirements for Safety
- **IEC 60601-1-2** – Top level standard for electromagnetic compatibility for electrical medical equipment
- **IEC 60601-2-38** – Requirements for safety of a hospital bed such as construction, mechanical strength, mechanical hazards, protection against fluid ingress, and essential performance
- **UL 60601** – Base standard for electromedical equipment
- **IPX 4** – Level of protection against ingress of fluid

For more information, contact your Hill-Rom representative or local distributor.

Your local Hill-Rom distributor is:

▶ Technical Data

Bed Length with Standard sleep deck	88.0" (223.5 cm)	Head Section Articulation0 through 72.5°
Bed Length with Special sleep deck	84.0" (213.5 cm)	Knee Section Articulation0 through 25°
Width (with siderails up)	43.7" (110 cm)	Caster Size	5.0" (12.7 cm)
Width (with siderails down)	36.6" (92 cm)	Safe Working Load	450 lbs. (205 kg)
Height of Sleep Deck (in low position)	18.3" (46.5 cm)	Operable Temperature Range*	+10° C to +40° C
Height of Sleep Deck (in high position)	30.1" (76.5 cm)	Operable Relative Humidity Range**	30% to 95% non-condensing
Underbed Clearance	5.5" (14.0 cm)		

* Product withstands storage/transportation in temperature range of -40° C to +70° C
 ** Product withstands storage/transportation in relative humidity range of 10% to 95%

Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2004 Hill-Rom Services, Inc.
 ALL RIGHTS RESERVED
 CTG674 04/30/04