


A HILL-ROM SOLUTION



HORIZON,
NEONATAL HORIZON,
UNIFLEX 1000, INTEGRIS 1000,
INTEGRIS 2000 HEADWALL SYSTEMS,
& EPICARE SUPPORT ARMS
ACCESSORIES

Hill-Rom
A HILLENBRAND INDUSTRY


INFUSION SUPPORT



P972D01 Utility/Infusion Mount, Dual Post

Height of post11"

Horizon[®]/Neonatal Horizon[®]/UniFlex[™]/Integris[™]



P971EF1 IV Hanger for Integris Patient Light



P158 Transfer Pole

Requires one P976H01 ISS Support Arm Mount.

Retracted position32"

Fully Extended57 1/2"

Diameter of top pole . . . 1 1/2"

Diameter of outer pole . . 1"

Horizon/Neonatal Horizon/UniFlex/Integris (Mounts to main raceway only on Integris)



P972N01 Short Pump Holder Pole

shown with P976 ISS Support Arm Mount.

Pole length24"

Pole diameter1/2"

Horizon/Neonatal Horizon/UniFlex/Integris




P984-00 Infusion Pump Holder for Horizon

with four hangers, requires two P102301 holders.

P984IN1 Infusion Pump Holder for Integris

with four hangers.



P971H00 IV Hanger

(shown with P1020A1 Utility/Bottle Slide)

Dimensions:

Outer2"W X 4 1/2"H

Inner1 1/4"W X 4"H

Horizon/Neonatal Horizon/EpiCare



P976H01 ISS Support Arm Mount for Horizon/Neonatal Horizon

Holds one P158 Transfer Pole.

P976IN1 ISS Support Arm Mount for UniFlex/Integris

Holds one P158 Transfer Pole.



P102302 Universal Holder


Horizon/Neonatal Horizon/UniFlex/EpiCare/Integris



P1020A1 Utility/Bottle Slide

Horizon/Neonatal Horizon/UniFlex/Integris/EpiCare

SHELVES AND BASKETS



P1014028 Utility Shelf


Laminate selection is required.

Dimensions:

Overall 20³/₈"W X 12¹/₈"D

Weight capacity 50 lbs.

Horizon/Neonatal Horizon




P1050EF1 Small Utility Shelf for Integris Patient Light

Laminate selection is required.

Dimensions:

Overall 12"W X 10"D

Weight capacity 15 lbs.




P105101 Large Storage Basket

Inner dimensions 16"W X 7¹/₂"H X 5"D

Functional loading up to 7 lbs.

Mounts to lower accessory track.

UniFlex/Integris 1000



P1029S3 Small Storage Basket-Stackable

Inner dimensions 21¹/₈"W X 4⁵/₁₆"H X 2⁷/₈"D

Functional loading up to 3 lbs.

Integris/EpiCare

P105002 Utility Shelf-Stackable

Laminate selection is required.

Dimensions:

Overall 20³/₈"W X 12¹/₈"D

Weight capacity 50 lbs.

UniFlex/Integris

P102102 Chart Holder-Stackable


Light Neutral Painted Steel

Dimensions 13"W X 8¹/₂"H X 2"D

Horizon/Neonatal Horizon/UniFlex/Integris

P1029S2 Small Storage Basket

Horizon/Neonatal Horizon/UniFlex



P1029L2 Large Storage Basket

Inner dimensions 16"W X 7¹/₂"H X 5"D

Functional loading up to 7 lbs.

Horizon/Neonatal Horizon/Integris/UniFlex/EpiCare

STORAGE



P101004B Wardrobe Right Hand/Wall Mount (shown)

P101005B Wardrobe Left Hand/Wall Mount

Patient Wardrobe with outer dimensions: 63 $\frac{1}{2}$ "H X 24"W at deepest end and 63 $\frac{1}{2}$ "H X 24"W X 9 $\frac{1}{8}$ "D

Laminate selection is required.

(Order for right hand hinged on left side of door, or left hand hinged on right side of door)

Horizon/Integris/UniFlex



P101107B Nursemate Door/Drawer-Horizon Mount

P101108B Nursemate Door/Drawer-Wall Mount

P1011IN1 Nursemate Door/Drawer-Integris Mount

Laminate selection is required.

Dimensions:
Inside drawer...
17"W X 3 $\frac{1}{4}$ "H X 13 $\frac{1}{8}$ "D
Overall...
21 $\frac{1}{4}$ "W X 26 $\frac{1}{4}$ "H X 16 $\frac{1}{4}$ "D
Overall weight...52 lbs.



P101137B Nursemate/Three Drawer-Horizon Mount

P101138B Nursemate/Three Drawer-Wall Mount

P1011IN2 Nursemate/Three Drawer-Integris Mount

Laminate selection is required.

Dimensions:
Inside drawer...
12 $\frac{1}{8}$ "W X 3 $\frac{1}{4}$ "H X 13 $\frac{1}{8}$ "D
Overall...
21 $\frac{1}{4}$ "W X 26 $\frac{1}{4}$ "H X 16 $\frac{1}{4}$ "D
Overall weight...52 lbs.




P101203B Mobile Storage Cart for Neonatal Horizon
Laminate selection is required.

Cushion is included.

Dimensions:
20"W X 23 $\frac{1}{8}$ "H X 17 $\frac{1}{8}$ "D

LIGHTS AND LAMPS



P963C01 Arm Lamp Incandescent Bulbs

P963C02 Arm Lamp Fluorescent or Incandescent Bulbs

P963M00 Arm Lamp Wall Mount

Arm extends to 46"
Requires power outlet and **P100801 Mount**

Horizon/Neonatal Horizon/UniFlex
Integris (not shown, mounts to main raceway only.)



P965-01 Exam Plus Light With Wall Mount
Extends 44". Requires power outlet and **P100901 Mount for Horizon/Neonatal Horizon**
P1009EF1 Mount for Integris Patient Light

Horizon/Neonatal Horizon/Integris

Also available:

PSPECLAX90 Exam Plus Light Ceiling Mount Kit with 15" rod. For ceiling heights from 8' to 8'6".

PSPECLAX91 Exam Plus Light Ceiling Mount Kit with 25" rod. For ceiling heights from 8'6"-9'6".



P102702 Utility Light

Horizon/Neonatal Horizon/
UniFlex/EpiCare/Integris

GASES



P1032A5 WYE Connector, Air
P1032AO5 WYE Connector, Oxygen
P1032A1 Port Cap, Air
P1032O1 Port Cap, Oxygen
P1032V1 Port Cap, VAC
P1032N1 Port Cap, Nitrous Oxide
Horizon/Neonatal Horizon/EpiCare



P1032V6 WYE Connector, VAC
Horizon/Neonatal Horizon/EpiCare



P2010 Integris Gas Block and Hose Assembly
Gas manufacturer and gas type (Quick or DISS) is required.




P1044 Horizon Gas Block and Hose Assembly
Gas manufacturer and gas type (Quick or DISS) is required.



P1045 Neonatal Horizon Gas Block and Hose Assembly
Gas manufacturer and gas type (Quick or DISS) is required.

MONITOR MANAGEMENT

See *Integris Monitor Arm brochure* for additional accessories



P100701 Horizon Monitor Slide (shown)

P1007N2 Neonatal Horizon Monitor Slide


P100702 Wall Mount Monitor Slide

P1007IN1 Integris Monitor Slide (needed only for P962-00)

P962-00 Tilt Swivel Monitor Arm (shown)

(When ordering, please specify monitor manufacturer and model #).

Horizon/Neonatal Horizon/Integris



P995P08 Satellite Rack & Mount-Hewlett Packard
Horizon/Integris/UniFlex

THERMOMETER MOUNT



P2027 Thermometer Mount-Diatek 6000/9000

UniFlex/Integris

P1020NB1X03 Thermometer Mount-Diatek 6000/9000

Horizon



P1020S1 Ophthalmoscope Mount (Extends Down)

Dimensions...
21¹/₁₆"W X 4¹/₄"H

Horizon/Neonatal
Horizon/UniFlex/Integris

OPHTHALMOSCOPE MOUNTS



P1020S2 Ophthalmoscope Mount, Welch Allyn 767

(Extends Down)
Dimensions...
6"W X 4¹/₄"H

Horizon/Neonatal
Horizon/UniFlex/Integris/
EpiCare




P1020S4 Ophthalmoscope Mount, Welch Allyn 767

(Extends Up)
Dimensions...
9¹/₁₆"W X 7"H X 2¹/₂"D

Horizon/Neonatal Horizon/
UniFlex/Integris/EpiCare

SPHYGMOMANOMETERS/CUFF BASKETS



P102503 Aneroid Sphygmomanometer Mount

P102802 Cuff Basket

P978-01 Aneroid Sphygmomanometer With Coil and Standard Size Cuff, Tycos

Horizon/Neonatal Horizon/UniFlex/Integris



P102402 Mercury Sphygmomanometer Mount, Baumanometer

P102802 Cuff Basket

Horizon/Neonatal Horizon/
UniFlex/Integris



P977-01 Mercury Sphygmomanometer With Coil and Standard Size Cuff, Baumanometer

P102403 Mercury Sphygmomanometer Mount

P102802 Cuff Basket


Neonatal Horizon/Integris

BUMPERS



P1016A2 Bed Bumper for Horizon

Dimensions 37¹/₄"W X 11¹/₄"H
 Extension from rail 3¹/₂"-1¹/₈"



P1016IN1 Bed Bumper for Integris

Dimensions 12"W X 15"H
 Extension from rail 3¹/₈"-1¹/₈"

MISCELLANEOUS ACCESSORIES



STAT Clock • Timer shown optionally flush-mounted in the Integris 2000 Headwall System. Must be specified when ordering for Integris/UniFlex. Must have provision in unit.

P101811 STAT Clock • Timer Pod Mount

Horizon/Neonatal Horizon
 Must have provision in unit.

P1018EF1 Integris STAT Clock • Timer Pod Mount

Integris Light Rail only.
 Must have provision in unit.

P967B00 STAT Clock • Timer, Panel Mount

Integris/UniFlex
 Must have provision in unit.

P967B01 STAT Clock • Timer, Wall Mount

Requires four gang masonry box 3¹/₂" deep provided by others. Includes pigtail for connection to building wiring.



**P103504 TV Arm Mount/
 RCA 1/4" Diameter
 for Horizon**

**P103505 TV Arm Mount/
 Bunting 1/2" Diameter
 for Horizon**

**P103506 TV Arm Mount/
 Phillips 1/4" Diameter
 for Horizon**



P103602 Cable Organizer
 Horizon/Neonatal Horizon/
 UniFlex/EpiCare/Integris



P2025EF1 Utility/Catheter Clip-Integris Patient Light
 Dimension 10"W



**P2026 Sharps Disposal
 Mount**
 Secures to rail via fastener.
 (When ordering please
 specify equipment
 manufacturer and model.)
 Horizon/UniFlex/Integris



**P101702 Waste
 Receptacle for Horizon**
 Dimensions...
 11¹/₈"W X 12¹/₈"H X 8¹/₈"D
**P1017IN1 Waste
 Receptacle for Integris
 2000 (not shown)**

MISCELLANEOUS ACCESSORIES



P1020N1 Extended Bottle Slide

Dimensions:

Outer2"W X 8¹/₈"H

Inner1¹/₄"W X 4"H

Horizon/Neonatal Horizon/
UniFlex/Integris



P1020B1 Bird Blender Mount (Extends Down)

Dimensions:

Outer . . .1²/₁₆"W X 4¹/₂"H

Inner2³/₄"W X 2¹/₈"H

Horizon/Neonatal Horizon/
UniFlex/Integris/EpiCare



P1020NB1 Bird Blender Mount for Neonatal Horizon (Extends Up)

Dimensions:

Outer . . .1²/₁₆"W X 4¹/₂"H

Inner2³/₄"W X 2¹/₈"H

Horizon/Neonatal Horizon



P1020HC1 Utility/Bottle Slide Mounts to lower accessory rail.

UniFlex/Integris 1000



P1020HC2 Utility/Catheter Clip

Mounts to lower accessory rail.

UniFlex/Integris 1000



P1020HC4 Utility Hook Mounts to lower accessory rail.

UniFlex/Integris 1000



P102602 Utility Hook With Mount

Horizon/Neonatal Horizon/
Integris/EpiCare



P102202 Resuscitator Bag Holder

Horizon/Neonatal Horizon/
UniFlex/Integris/EpiCare

DIRECT CUSTOMER SERVICE HOTLINE:
USA 800-445-3730

Hill-Rom[®]

A HILLERBRAND INDUSTRY

Batesville, IN 47006 USA • FAX 812-934-1697
International • 812-934-8173 • FAX 812-934-7191
www.hill-rom.com

CTG036ra

Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the expressed written warranty extended on the sale or rental of its products. ©1999 Hill-Rom Company, Inc. ALL RIGHTS RESERVED.