



## Wheelchair

- Anti-tip mechanism helps prevent the chair from tipping backwards
- Adjustable width leg supports to better accommodate patient size
- Attached IV pole and rear push bar let caregivers transport patients more efficiently
- Full-length removable, height-adjustable arm rests
- Metal mag-style 24" rear wheels and semi-pneumatic 8" x 2" front caster wheels help provide dependable support for mobility
- Dimensions:

<i>Height</i>	Seat 18.75" Overall 37"
<i>Width</i>	Seat 30" Overall 39"
<i>Depth</i>	Seat 20" Overall 33"
<i>Weight</i>	Net Weight 88 pounds Maximum recommended patient weight 800 pounds



## Walker

- Adjustable height to accommodate patient height
- Anodized double lumen aluminum tubing for increased strength
- Double braces on frame provide additional support
- Foam hand grips help to provide a secure surface
- Dimensions:

<i>Height</i>	Adjusts 30" to 40"
<i>Width</i>	Top 25.5" Base 27"
<i>Depth</i>	Top 14.5" Base 19"
<i>Weight</i>	Net Weight 10 pounds Maximum recommended patient weight 800 pounds



## Shower/Commode Chair

- Rubber stops on all legs help to secure placement
- Commode lid acts as a backrest for enhanced patient comfort
- Adjustable height allows for placement over hospital wall mounted toilets
- Stainless steel construction allows for use in the shower
- Dimensions:

<i>Height</i>	Seat Adjusts 18" to 24" Overall 27.75" to 33.75"
<i>Width</i>	Seat 28" Overall 30"
<i>Depth</i>	Seat 22.5"
<i>Weight</i>	Net Weight 54 pounds Maximum recommended patient weight 800 pounds



## Body Mass Index Table

To use the table, find the appropriate height in the left-hand column labeled Height. Move across to a given weight. The number at the top of the column is the BMI at that height and weight. Pounds have been rounded off.

### Extreme Obesity

BMI	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54
Ht. (in)	Body Weight (pounds)																		
58	172	177	181	186	191	196	201	205	210	215	220	224	229	234	239	244	248	253	258
59	178	183	188	193	198	203	208	212	217	222	227	232	237	242	247	252	257	262	267
60	184	189	194	199	204	209	215	220	225	230	235	240	245	250	255	261	266	271	276
61	190	195	201	206	211	217	222	227	232	238	243	248	254	259	264	269	275	280	285
62	196	202	207	213	218	224	229	235	240	246	251	256	262	267	273	278	284	289	295
63	203	208	214	220	225	231	237	242	248	254	259	265	270	278	282	287	293	299	304
64	209	215	221	227	232	238	244	250	256	262	267	273	279	285	291	296	302	308	314
65	216	222	228	234	240	246	252	258	264	270	276	282	288	294	300	306	312	318	324
66	223	229	235	241	247	253	260	266	272	278	284	291	297	303	309	315	322	328	334
67	230	236	242	249	255	261	268	274	280	287	293	299	306	312	319	325	331	338	344
68	236	243	249	256	262	269	276	282	289	295	302	308	315	322	328	335	341	348	354
69	243	250	257	263	270	277	284	291	297	304	311	318	324	331	338	345	351	358	365
70	250	257	264	271	278	285	292	299	306	313	320	327	334	341	348	355	362	369	376
71	257	265	272	279	286	293	301	308	315	322	329	338	343	351	358	365	372	379	386
72	265	272	279	287	294	302	309	316	324	331	338	346	353	361	368	375	383	390	397
73	272	280	288	295	302	310	318	325	333	340	348	355	363	371	378	386	393	401	408
74	280	287	295	303	311	319	326	334	342	350	358	365	373	381	389	396	404	412	420
75	287	295	303	311	319	327	335	343	351	359	367	375	383	391	399	407	415	423	431
76	295	304	312	320	328	336	344	353	361	369	377	385	394	402	410	418	426	435	443

Hill-Rom is a registered trademark of Hill-Rom Services, Inc.  
Hill-Rom reserves the right to make changes without notice  
in design, specifications and models. The only warranty  
Hill-Rom makes is the express written warranty extended on  
the sale or rental of its products.

© 2003 Hill-Rom Services, Inc.  
ALL RIGHTS RESERVED.

CTG128 Rev 2 04/21/05