

Bedside Cabinets - Classic

PRODUCT SUMMARY



The Bedside Cabinet provides ease of storage where needed the most.

CLASSIC 3-DRAWER BEDSIDE CABINET
(BSC-100-r3o-rok)



Ultra-tough tops are stain and scratch resistant.



Easy-roller casters allow unit to be moved for patient and staff convenience.



Generously rounded corners of storage drawers make cleaning easier.

CLASSIC BEDSIDE CABINET



Replaceable
Panels

Ventilated
Back

Steel Frame

Steel Base



SPECIFICATIONS:

Overall Height	31.5" / 80 cm
Overall Width	21.5" / 55 cm
Overall Depth	17.5" / 44 cm
Shipping Weight	66 lbs / 30 kg
Floor Print	21.5" W x 17.5" D (55 cm x 44 cm)

Hill-Rom
A HILLENBRAND INDUSTRY

Batesville, IN 47006 USA • 812-934-7777 • Fax 812-934-8189
Canada • 905-206-1355 • Fax 905-206-0561
International • 812-934-8173 • Fax 812-934-7191
www.hill-rom.com

CTG159ra

Due to limitations of the printing process, colors may vary slightly. When specifying, obtain color samples by contacting the fabric or finish manufacturer directly.

Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
©2002 Hill-Rom Services, Inc.
ALL RIGHTS RESERVED.