

Bedside Cabinets - Elite

PRODUCT SUMMARY



The Bedside Cabinet provides ease of storage where needed the most.

ELITE 3 - DRAWER BEDSIDE CABINET
(BSC300-r30-rok)



Galley rails help prevent patient's items from sliding off the top surface.



Pullout tray offers extra work space for caregivers.



Easy-roller casters allow unit to be moved for patient and staff convenience.



Generously rounded corners of storage drawers make cleaning easier.

Hill-Rom
A HILLENBRAND INDUSTRY

ELITE BEDSIDE CABINET



Replaceable Panels

Ventilated Back

Steel Frame

Steel Base



SPECIFICATIONS:

Overall Height	31.5" / 80 cm
Overall Width	21.5" / 55 cm
Overall Depth	17.5" / 44 cm
Shipping Weight	80 lbs / 37 kg
Floor Print	21.5" W x 17.5" D (55 cm x 44 cm)

Hill-Rom
A HILLENBRAND INDUSTRY

Batesville, IN 47006 USA • 812-934-7777 • Fax 812-934-8189
Canada • 905-206-1355 • Fax 905-206-0561
International • 812-934-8173 • Fax 812-934-7191
www.hill-rom.com

CTG161ra

Due to limitations of the printing process, colors may vary slightly. When specifying, obtain color samples by contacting the fabric or finish manufacturer directly.

Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
©2002 Hill-Rom Services, Inc.
ALL RIGHTS RESERVED.