

Hill-Rom Power Columns



Power Column

At a glance

- Four-sided access to the patient
- Efficient use of premium space
- Two-sided accessory tracks on all four corners of the column
- Numerous accessories available
- Two sizes of power columns (23" and 15") for maximum flexibility

Maximum access, maximum flexibility

Hill-Rom power columns offer the best of both worlds: 360° access to critically ill patients, and the best use of space in the crowded acute care environment.

Hill-Rom power columns allow caregivers to move freely around the bed, including the critical area around the head of the patient. Free of the wall, the bed can be pivoted around to allow the patient to see out the window or watch television, without disturbing vital services.


Hill-Rom Power Columns



Built and configured just for you

Each power column is designed and built to individual specifications. The unique accessory track system provides full-length, two-sided tracks on all four corners, allowing staff to position equipment where it's needed most.

Both power columns offer full service capacity, making it easy for staff to focus on patient care.



Small Power Column
Side View



Small Power Column
Front View

Technical data

Column cross-section

Power Column 23 1/8" W x 7" D
 Small Power Column 15 1/8" W x 7" D

Floor Space (Base)

Power Column 21 1/2" W x 6" D
 Small Power Column 13 1/2" W x 6" D

Ceiling Space

Power Column 23 1/4" W x 10 1/2" D
 Small Power Column 15 1/4" W x 10 1/2" D

Height (Floor to Ceiling) 8'0"–9'6" in 1" increments

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2003 Hill-Rom Services, Inc
 ALL RIGHTS RESERVED.

CTG267rb 10/21/03

