

Affinity® Three Birthing Bed Accessories





P3613 LABOR BAR
Allows the mother to assume a squatting position.
Load capacity 100lbs (45kg)



P451 FOOT REST BAR
Provides a convenient foot rest for the physician
Load capacity 150lbs (68kg)



P3704 CALF SUPPORTS
Mounted to back of foot supports for easy access.
Load capacity 60lbs (27kg)



P35745A TELESCOPING CALF SUPPORTS
Fully-padded with telescoping height adjustment.
Load capacity 60lbs (27kg)



P7625C LEG SUPPORTS
Padded full leg support with ball joint and pulling handles
Load capacity 60lbs (27kg)



P7634C TELESCOPING LEG SUPPORTS
Padded full leg supports with with telescoping height adjustment.
Load capacity 60lbs (27kg)



P3616C CALF SUPPORTS

Designed specifically for the Affinity® I and II beds. Comes standard with plastic holster for permanent attachment.



P3624E SIDERAIL PADS

The pads, designed to fit Affinity® I and II birthing beds, provide a soft cover for the rails while maintaining access to patient nurse controls.

AFFINITY® BIRTHING BED ACCESSORIES



HEAD SECTION SLIP COVERS

- P3625 V-Cut, without air
- P3626 Straight Edge, without air
- P3627 Straight Edge, with air
- P3628 V-Cut, with air

FOOT SECTION SLIP COVERS

- P3630 V-cut
- P3631 Straight Edge

Exclusive form-fitting Dartex® Slipcover, with water proof seams, extends the durability of the mattress surface and protects inner foam from tears, punctures or general leakage. The Ultra-Fresh anti-microbial, latex-free Dartex® cover is made from a stretchable polyurethane which is soft, durable and easy-to-clean.

AFFINITY® BIRTHING BED HEADBOARD OPTIONS




P3629 BLOW-MOLDED HEADBOARD WITH LAMINATE




P3617 CONTEMPORARY WOOD HEADBOARD




P93-3617 FREEDOMHILL™ COLLECTION WOOD HEADBOARD




P92-3617 LIBERTYHILL™ COLLECTION WOOD HEADBOARD



P3618 PENDANT
P497 PENDANT HOLDER-FOR SIDERAIL
P5362 PENDANT KIT (P3618, P497)



P3619 WEDGE ASSEMBLY
Converts V-Cut to straight edge quickly and easily to accommodate physician preference.



P3623 FOLEY HOOK KIT
Load capacity 20lbs (9kg)



P9625 ANESTHESIA SCREEN
Load capacity 5lbs (2.3kg)



P9620 ANESTHESIA ARM BOARD
Length: 39.25"(99.7cm), Height: 6"(15.24cm),
Thickness: .25"(.64cm)
Load capacity 30lbs (14kg)




P278 INSTRUMENT TRAY
Dimension, outside: 22.63" (57.5cm) x
9.25" (23.5cm)
Load capacity 25bs (11kg)



P495 PATIENT PHONE
P496 PATIENT PHONE ADAPTER-
SIDERAIL



P2221-01 PERMANENT IV POLE
Load capacity: 25 lbs (11kg)
Mounts to the head section of the bed. Can be
used on either side.



P157-10 COMPLETE INFUSION SUPPORT SYSTEM
P159 OFFSET BAR
P158 TRANSFER POLE
Load capacity: 40 lbs (18kg)



P156 ISS MOBILE STAND SYSTEM
P156P MOBILE STAND TRANSFER POLE
P156B MOBILE STAND BASE WITH CASTERS




P2217 STANDARD IV POLE
Load capacity: 25 lbs (11kg)



P3706 SIDERAIL PADS FOR AFFINITY® THREE BIRTHING BED
Soft cover for the rails while maintaining access to patient nurse controls.



P3612 COMFORT PAD
Designed with 2" (5.1cm) egg crate foam to provide an extra layer of comfort.



B7831 DELIVERY DRAPES
Case of 20 Full Size 32" x 58" (81.28cm x 147.3cm)
SA6603 DELIVERY DRAPES
Case of 40 Half Size 32" x 29" (81.28cm x 73cm)



P276 O₂ TANK HOLDER
Holds standard O₂ tank.
Load capacity 30lbs (14kg)



SA4585 SHEET SET-V-CUT SURFACE
SA4586 SHEET SET-STRAIGHT SURFACE
1528 BEDDING SET-V-CUT SURFACE
1529 BEDDING SET-STRAIGHT SURFACE

ADDITIONAL PRODUCTS



P118B524 SURGICAL STOOL
Adjusts from 20" to 24" (50.8cm to 61cm) for convenience and comfort using a foot controlled gas spring.
15" Seat Diameter (38.1cm)
Load capacity 300lbs (136kg)



P7925 PRIMA PROCEDURAL LIGHT
Designed to provide no-drift, shadow-free surgically correct lighting. Mounted on a four-leg base with two locking casters for stability. Base 23"D x 19"W x 60"-73"H (58.4cm x 48.3cm x 152.4-185.4cm)



P15-05 BIRTH MIRROR
Adjustable so patient can view the birth from almost any delivery position. 24"W x 62½"H (609.6cm x 158.8 cm)
Mirror 14"W x 18"H (35.6cm x 45.7cm)

ACCESSORIES

Anesthesia Screen	P9625
Arm Board	P9620
Calf Supports, Affinity I, II	P3616C
Calf Supports, Affinity Three.....	P3704
Calf Supports-Telescoping	P35745A
Comfort Pad	P3612
Disposable Drapes-Full Size	B7831
Disposable Drapes-Half Size	SA6603
Foley Hook Kit	P3623
Foot Rest Bar.....	P451
Full Leg Supports	P7625C
Full Leg Supports Telescoping	P7634
Instrument Tray	P278
Complete ISS System	P157-10
Offset Bar	P159
Transfer Pole	P158
Mobile Stand Transfer Pole.....	P156P
Mobile Stand Base with Casters	P156B
Standard IV Pole.....	P2217

Labor Bar	P3613
O ₂ Tank Holder.....	P276
Pendant Control Holder-Siderail	P497
Permanent IV Pole	P2221-01
Patient Phone	P495
Patient Phone Adaptor	P496
Sheet set-V-cut Surface.....	SA4585
Sheet set-Straight Surface	SA4586
Bedding set-V-cut Surface	SA1528
Bedding set-Straight Surface	SA1529
Slip Cover-Head Section	P3625
	P3626
	P3627
	P3628
Slip Cover-Foot Section Straight	P3631
V-Cut.....	P3630
Siderail Pads, Affinity I & II.....	P3624
Siderail Pads, Affinity Three	P3606
Wedge Assembly	P3619

HOSPITAL INSTALLED OPTIONS

- 8" (20 cm) Casters
- Air Support System
- SideCom® Communication Controls
- Variety of Wood Headboards
- Foot Section Conversion Kit
- Seat Section Conversion Kit

Comfort Pads and siderail pad approved for use in San Francisco (CAL 129 accepted in lieu of CAL 121), Maryland (CAL 129), Boston (IX-11), North Carolina, Ohio and Washington, DC (with Dartex)

DIRECT CUSTOMER SERVICE HOTLINE:
USA 800-445-3730
CANADA 800-267-2337

Hill-Rom

A HILLENBRAND INDUSTRY

Batesville, IN 47006 USA • 812-934-7777 • Fax 812-934-8189
Canada • 905-206-1355 • Fax 905-206-0561
International • 812-934-8173 • Fax 812-934-7191
www.hill-rom.com

CTG291

Dartex® is a registered trademark of Courtaulds Textiles (Holdings) Limited.

Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
©2001 Hill-Rom Services, Inc.
ALL RIGHTS RESERVED.