

The Impression™ Series Headwall Collection





Today's competitive marketplace requires increased emphasis on aesthetics throughout all areas of the patient environment. The Impression™ Series Headwalls from Hill-Rom, offer a solution to the challenge of creating a warm relaxing atmosphere for patients and their visitors without compromising the scope or placement of service.

An indirect wall light provides soft lighting to enhance the room's overall appearance and create a warmer atmosphere for patients and visitors. An integrated light provides abundant lighting directed in front of the patient for reading.

The recessed service side of the cabinet is deep enough to provide ample space for most secondary equipment such as suction canisters and tubing. Equipment may be set up and ready for use behind the closed tambour doors, thus saving valuable set-up time when needed in emergency situations. Vacuum bottle slide holders are stored out of sight, yet easily accessible when needed.

The headwall also features a safety interlock switch that prevents the bed from further movement when contact is made with the light or the lower frame of the unit.

While beautifully enhancing the room decor, the Impression™ Series Headwall delivers all of the services that moderate-acuity patients require. And, room furnishings from Hill-Rom may be fully coordinated to optimize aesthetic appeal.

Of course, Impression™ Series Headwalls are just one solution from Hill-Rom that meets health care environment needs around the world.

The exterior trim is made with solid oak^(R) or maple^(M) and comes in your choice of finishes:



Medium Oak - 82R



Wild Cherry - T9M



Royal Oak
H5M or H5R



Dark Cherry - H4M



Bleached Oak - J2R



Light Oak - R3R



Natural
31R or 31M

Due to limitations of the printing process, colors may vary slightly. When specifying, obtain color samples by calling our direct customer service hotline.

Baby Impression 783D

This unit is ideal for low acuity areas such as med surg, outpatient, and LDR/P. It is designed to save space, yet houses necessary services.

Matching units without services are also available to provide additional storage of nursing supplies while maintaining the complete, coordinated aesthetics of the room.



Concealed Services

Tambour doors easily slide open to reveal such necessary services as electric, gases, nurse call, sphygmomanometer and IV holders. This unit is shown with DiamondCare[®] outlets, from Hill-Rom, which are available in all popular key styles.

Impression™ Series Dimensions

	PRIVATE MODEL 780D	SEMI-PRIVATE MODEL 780D	BABY IMPRESSION MODEL 783D
Overall			
Width	7'10¼" (2,39 m)	14'½" - 14'9½" (4,28 m - 4,51 m)	37" (94 cm)
Height	26" (66,0 cm)	26" (66,0 cm)	26" (66,0 cm)
Depth			
<i>Of Service Wings</i>	12" (30,5 cm)	12" (30,5 cm)	12" (30,5 cm)
<i>At Center Access Panel</i>	5¼" (13,3 cm)	5¼" (13,3 cm)	
Weight	215 lbs. (97,5 kg)	410 lbs. (186 kg)	64 lbs. (29 kg)
Inside Service Wing			
Depth of Plain Type for Full Height	6¼" (15,9 cm)	6¼" (15,9 cm)	6¼" (15,9 cm)
Depth of Recess Type for Full Height	9¼" (4,1 cm)	9¼" (4,1 cm) (extending up to 11¼") (28,6 cm)	N/A

DIRECT CUSTOMER SERVICE HOTLINE:
USA 800-445-3730
CANADA 800-267-2337

Hill-Rom
A HILLENBRAND INDUSTRY

Batesville, IN 47006 USA • 812-934-7777 • Fax 812-934-7038
Canada • 905-206-1355 • Fax 905-206-0561
International • 812-934-8173 • FAX 812-934-7191
www.hill-rom.com

CTG298



Healthcare Design Research Alliance™ Mission Statement

Provide a forum for committed manufacturers and selected others to collaborate in the research, design and marketing of innovative products and services that improve the quality of healthcare interiors and contribute to a leadership image and the business objectives of each HDRA participant.



Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
©2001 Hill-Rom Services, Inc.
ALL RIGHTS RESERVED.