

The ImageMakers[®]
Headwall Collection





All furniture pieces and headwalls shown are finished in Natural Maple.

Hill-Rom offers a beautiful solution to the problem of how to create a warm, relaxing atmosphere for patients without compromising the scope or placement of services.

The ImageMakers® Headwall Collection delivers all of the services to moderately acute patients while greatly enhancing a room's decor. Moreover, the depth of the service cabinets allows most secondary equipment to be attached and stored behind closed doors, hidden but ready for use.

Each headwall features deep polished wood, hand-rubbed for clarity and a high luster sheen. The finish is resistant to water, alcohol and nail polish remover.

Both the full size and small size headwalls are available in three distinct styles and a variety of decorator stains. And, each will integrate beautifully with Hill-Rom's ImageMakers™ and Perinatal Furniture Collections to optimize aesthetic appeal and create a coordinated systems approach to your patient room.

Available in your choice of finishes:



Medium Oak - 82R



Wild Cherry - T9C



Royal Oak - H5R



Dark Cherry - H4C



Bleached Oak - J2R



Light Oak - R3R

Due to limitations of the printing process, colors may vary slightly. When specifying, obtain color samples by calling our direct



The Contemporary Headwall

Features solid oak and veneers crafted to take advantage of the beauty and utility of the best in modern design.



The FreedomHill™ Headwall

Features solid oak and veneers to create an environment of distinction and timeless beauty.



The LibertyHill™ Headwall

Features solid cherry and veneers crafted to offer the old-world charm of 18th Century Traditional furniture.



Storage Unit

To maintain a coordinated look to the room, a matching storage unit for your nursing supplies is available in the same styles and stains as the headwall unit.



Small Size Headwall

The Small Size Headwall is ideally suited for low acuity areas such as Sub-Acute Care, Out Patient, Rehab, etc. and is available in the same styles and stains as the full size unit.

ImageMakers® Dimensions

	FULL SIZE • MODEL 790	SMALL SIZE • MODEL 793
Overall		
Width	96" (243.84 cm)	28½" (72.39 cm)
Height	26½" (66.36 cm)	26½" (66.36 cm)
Depth	14½" (35.88 cm)	14½" (35.88 cm)
Weight	190lbs. (86 kg)	79lbs. (36 kg)
Inside Cabinet		
Width	24" (60.96 cm)	24" (60.96 cm)
Height	21¾" (55.25 cm)	21¾" (55.25 cm)
Depth - to device panel		
<i>Doors Open</i>	8¼" (20.96 cm)	8¼" (20.96 cm)
<i>Doors Closed</i>	7½" (19.05 cm)	7½" (19.05 cm)
Depth - to back wall		
<i>Doors Open</i>	11¾" (29.85 cm)	11¾" (29.85 cm)
<i>Doors Closed</i>	11" (27.94 cm)	11" (27.94 cm)

NOTE: The ImageMakers® Headwall is also available in models designed to cover in-wall services.

DIRECT CUSTOMER SERVICE HOTLINE:

USA 800-445-3730

CANADA 800-267-2337

Hill-Rom®

A HILLENBRAND INDUSTRY
 Batesville, IN 47006 USA • 812-934-7777 • Fax 812-934-7038
 Canada • 905-206-1355 • Fax 905-206-0561
 International • 812-934-8173 • FAX 812-934-7191
 www.hill-rom.com

CTG299



Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
 ©2001 Hill-Rom Services, Inc.
 ALL RIGHTS RESERVED.