

Uniflex[®] 1000 Headwall System



Uniflex® 1000 Headwall System

A Practical Solution to Equipment and Service Management

The UniFlex® 1000 Headwall, from Hill-Rom, brings a clean, contemporary look to intermediate care areas of the hospital.

This efficient headwall provides an alternative to in-wall construction systems. It is delivered pre-plumbed and pre-wired to permit single point connections for each service type and mounts directly to the wall with brackets that are provided.

Available with an optional vertical chase, the UniFlex® headwall organizes vacuum, oxygen, air, electrical and communications within easy reach of the staff. Ancillary equipment is easily attached along the upper and lower accessory tracks allowing the system to adapt to changing care requirements.

The headwall comes in variable lengths and a wide choice of laminate colors. It may be used alone, or as a system with the addition of the Specialite™ and Bed Locator from Hill-Rom.

When an economical solution is demanded, the UniFlex® 1000 Headwall is the answer.



Dimensions

Length (1 rail) . . .3' minimum up to 9' maximum
(91.44 cm up to 274.32 cm)

For units longer than 9',
rails easily connect together

Rail Height14" (35.56 cm)

Rail Depth3¹⁵/₁₆" (10 cm)

Rail Distance from floor to the device
centerline48" to 60" (standard)
(121.92 cm x 152.4 cm)

Vertical Chase . . .5⁵/₈" deep by 13³/₄" wide
(14.29 cm x 34.93 cm)

Hill-Rom®

A HILLENBRAND INDUSTRY

Batesville, IN 47006 USA • 812-934-7777 • Fax 812-934-7038
Canada • 905-206-1355 • Fax 905-206-0561
International • 812-934-8173 • Fax 812-934-7191

www.hill-rom.com

CTG309

Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
©2001 Hill-Rom Services, Inc.
ALL RIGHTS RESERVED.