

# TranStar®

## *Procedural Stretcher*



Steering Plus™ System\* for superior maneuverability and sharp cornering

One-Step® TuckAway Siderails\* for zero-gap transfers

Ergonomic Push Handles\*\* to reduce caregiver strain during transport

*\*Standard on PC-300 through PC-550 model procedural stretchers*

*\*\*Standard on PC-200 through PC-550 model procedural stretchers*

**Hill-Rom®**  
A HILLENBRAND INDUSTRY

# TranStar® Procedural Stretcher

The TranStar® Procedural Stretcher delivers new levels of mobility, comfort and efficiency. Designed with the patient and staff in mind, the unit offers such features as a unique steering system, custom siderails, ergonomic transport handles, a low Lo height, radiolucent surface and convenient, functional controls.

These features, combined with its full array of options and accessories, make the TranStar® Procedural Stretcher the ideal solution for patient care and transport throughout the health care facility.

## Procedural Stretchers Value Package Configurations

Model	26" Width	30" Width	Ergonomic Push Handles	IV Pole(s)	5th Wheel Steering	Knee Gatch	Radio-lucent Surface	Footend Pedals	Mattress Options				
									3" Thick	4" Thick	5" Thick	26" Width	30" Width
PC-550		x	x	2	x	x	x	x	Std	Op	Op		Std
PC-500	x		x	2	x	x	x	x	Std	Op	Op	Std	
PC-450		x	x	1	x	x			Std	Op	Op		Std
PC-400	x		x	1	x	x			Std	Op	Op	Std	
PC-350		x	x	1	x				Std	Op	Op		Std
PC-300	x		x	1	x				Std	Op	Op	Std	
PC-200	x		x	1					Std	Op	Op	Std	
PC-100	x								Std	Op	Op	Std	

Note: KX1 Added Coverage Siderail, with footend gap, is available for an additional cost on model PC-300 through PC-550 stretchers.

CAL 129 Flammability requirement is an option on all mattresses.

Labeling: **Standard decals include:** Emergency, Endoscopy, Cath Lab, GI Lab, PACU, Transport, Maternity, Radiology, Surgery, and Ambulatory Surgery.

**Custom decals:** See your Hill-Rom Customer Service Agent or Account/Territory Manager if custom decals are desired.

### Technical Data PC-100–PC-200

Low Position — Floor to Top of Surface: .....	23" (58.4 cm)
High Position — Floor to Top of Surface: .....	35" (88.9 cm)
Overall Length: .....	83" (210.8 cm)
Overall Width: .....	32" (81.3 cm)
Siderail Length: .....	47" (119.4 cm)
Siderail Height Above Mattress/Surface: .....	14" (36.8 cm)
Mattress Size: .....	26" (66 cm) x 75" (190.5 cm)
Maximum Head Elevation: .....	90°
Maximum Trendelenburg/Reverse Trendelenburg: .....	18°
Floor to Base Clearance: .....	3.5" (8.9 cm)
Caster Size: .....	8" (20.3 cm)
Maximum Weight Limit: .....	500 lbs. (226.8 kg)

### Technical Data PC-300–PC-550

Low Position — Floor to Top of Surface: .....	21.5" (54.6 cm)
High Position — Floor to Top of Surface: .....	35" (88.9 cm)
Overall Length: .....	83" (210.8 cm)
Overall Width — Siderails Up: .....	32" (81.3 cm) or 36" (91.4 cm)
Overall Width — Siderails Stored: .....	26.5" (67.3 cm) or 30.25" (76.8 cm)
Siderail Length: .....	47" (119.4 cm)
Siderail Height Above Mattress/Surface: .....	14.5" (36.8 cm)
Mattress Size: .....	26" (66 cm) x 75" (190.5 cm)
Maximum Head Elevation — Manual Inclination: .....	90°
Maximum Knee Elevation: .....	25°
Maximum Trendelenburg/Reverse Trendelenburg: .....	18°
Floor to Base Clearance: .....	3.5" (8.9 cm)
Caster Size: .....	8" (20.3 cm)
Maximum Weight Limit: .....	500 lbs. (226.8 kg)

Direct Customer Service Hotline: USA 800-445-3730 | Canada 800-267-2337

**Hill-Rom**  
A HILLENBRAND INDUSTRY

Batesville, IN 47006 USA • 812-934-7777 • Fax 812-934-8189  
Canada • 905-206-1355 • Fax 905-206-0561  
International • 812-934-8173 • Fax 812-934-7191

www.hill-rom.com

CTG465ra



Patent Pending

PC100 & PC200 Models Are Not CE Approved

Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2002 Hill-Rom Services, Inc.  
ALL RIGHTS RESERVED.

## Features

- Ideal for use in Emergency, OR, PACU, and Outpatient
- Low height improves on/off transfers and increases patient safety
- Durable Siderails quickly tuck away for zero-gap transfers and full caregiver access (PC-300–PC-550 models)
- Ergonomically correct Transport Handles reduce caregiver strain during patient transport (PC-200–PC-550 models)
- TRI-PEDAL Foot Controls allow even lowering of patients
- TranStar® Stretchers feature a shock-absorbing frame and siderail system for improved patient comfort and satisfaction.

## Accessories

- IV Transporter
- Patient Tray
- Utility Tray
- Restraint Straps
- Oxygen Tank Holder
- Liquid Oxygen Tank Holder
- Infusion Support System
- Siderail Covers
- Armboard
- Detachable Footboard
- Mattress Pad
- Lateral X-ray Cassette Holder
- Chart Holder
- Convertible Footboard