

TranStar[®]

OB/GYN Stretcher



Integrated Foot Supports save time by eliminating the need to search for stirrups or calf supports

Sliding litter, with store-away foot section, enables quick and easy patient positioning

Steering Plus... System for superior maneuverability and sharp cornering

OneStep[®] TuckAway Siderails for zero-gap transfers

Ergonomic Push Handles to reduce caregiver strain during transport

Hill-Rom[®]
A HILLENBRAND INDUSTRY

TranStar®

OB/GYN Stretcher

Hill-Rom's TranStar® OB/GYN Stretcher combines innovation, mobility, and ease of use into a lightweight and highly durable unit. Designed for use in Maternity and Emergency, the OB/GYN Stretcher offers an array of features found only on the TranStar® line of Stretchers.

The Stretcher features integrated foot supports that eliminate searching for stirrups. The convenient store-away foot section eliminates heavy lifting and searching for storage space. And, the Stretcher's litter slides toward the foot end, eliminating the need to reposition the patient. The OB-550 package provides an integrated fiber optic exam light designed for precision and convenience. The TranStar® OB/GYN Stretcher is the ideal solution for pelvic examinations.

OB/GYN Stretchers Value Package Configurations

Model	30" Width	Ergonomic Push Handles	IV Pole(s)	5th Wheel Steering	Radio-lucent Surface	Factory Mounted Light	Mattress Options	
							3" Thick	30" Width
OB-550 OB/GYN	x	x	1	x	x	x	Std	Std
OB-450 OB/GYN	x	x	1	x	x		Std	Std

Labeling: **Standard decals include:** Emergency, Endoscopy, Cath Lab, GI Lab, PACU, Transport, Maternity, Radiology, Surgery, and Ambulatory Surgery

Custom decals: See your Hill-Rom Customer Service Agent or Account/Territory Manager if custom decals are desired.

Technical Data

Low Position — Floor to Top of Surface:	24.25" (61.6 cm)
High Position — Floor to Top of Surface:	38" (96.5 cm)
Overall Length:	80.5" (204.5 cm)
Overall Width — Siderails Up:	35.5" (90.2 cm)
Overall Width — Siderails Stored:	31" (78.7 cm)
Siderail Length:	37" (94.0 cm)
Siderail Height Above Surface:	13" (33.0 cm)
Mattress Size:	30" x 76.5" (76.2 x 194.3 cm)
Maximum Head Elevation:	90°
Maximum Trendelenburg/Reverse Trendelenburg:	18°
Floor to Base Clearance:	3.5" (8.9 cm)
Caster Size:	8" (20.3 cm)
Maximum Weight Limit:	500 lbs. (226.8 kg)

Direct Customer Service Hotline: USA 800-445-3730 | Canada 800-267-2337

Hill-Rom®

A HILLENBRAND INDUSTRY

Batesville, IN 47006 USA • 812-934-7777 • Fax 812-934-8189

Canada • 905-206-1355 • Fax 905-206-0561

International • 812-934-8173 • Fax 812-934-7191

www.hill-rom.com

CTG466



Patent Pending

Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2002 Hill-Rom Services, Inc.
ALL RIGHTS RESERVED.

Features

- Ideal for use in Emergency and Maternity
- Store-Away Foot Section eliminates heavy lifting and storage concerns
- Sliding Litter enables quick and easy patient positioning
- Integrated Foot Supports provide unmatched convenience and save time by eliminating the need to search for stirrups and calf supports
- Foot Supports telescope upward and outward in multiple positions. Attached straps and built-in heel supports provide secure positioning of the feet
- Integrated Catch Basin, combined with patient grip handles located on the siderails, offers convenience and saves time
- TranStar® Stretchers feature the GentleRide™ shock-absorbing frame and siderail system for improved patient comfort and satisfaction.

Accessories

- IV Transporter
- Patient Tray
- Restraint Straps
- Oxygen Tank Holder
- Liquid Oxygen Tank Holder
- Infusion Support System
- Siderail Covers
- Armboard
- Headboard
- Upright X-ray Cassette Holder
- Integrated Telescoping Calf Supports
- Catch Basin