

# TranStar<sup>®</sup>

## *Trauma Stretcher*



Full-length, head-to-toe, X-ray cassette positioning capability  
saves time and reduces patient risk

Steering Plus<sup>™</sup> System for superior maneuverability and sharp cornering

OneStep<sup>®</sup> TuckAway Siderails for zero-gap transfers

Ergonomic Push Handles to reduce caregiver strain during transport

**Hill-Rom**<sup>®</sup>  
A HILLENBRAND INDUSTRY

# TranStar® Trauma Stretcher

The Hill-Rom TranStar® Trauma Stretcher is a highly mobile and durable unit that is ideal for trauma situations when performance is paramount.

The unit features full-length X-ray positioning with an upright chest X-ray cassette holder. The cassettes can be positioned from head-to-foot and can be raised to the sleep surface for optimal image quality.

In addition, an optional lateral X-ray accessory is available.

## Trauma Stretcher Value Package Configuration

Model	26" Width	Ergonomic Push Handles	IV Pole(s)	5th Wheel Steering	Radio-lucent Surface	Footend Pedals	Mattress Options			
							3" Thick	4" Thick	5" Thick	26" Width
TR-500 Trauma	x	x	1	x	x	x	Std	Option	Option	Std

Labeling: **Standard decals include:** Emergency, Endoscopy, Cath Lab, GI Lab, PACU, Transport, Maternity, Radiology, Surgery, and Ambulatory Surgery

**Custom decals:** See your Hill-Rom Customer Service Agent or Account/Territory Manager if custom decals are desired.

## Technical Data

Low Position — Floor to Top of Surface: .....	25.62" (65.1 cm)
High Position — Floor to Top of Surface: .....	38.5" (97.8 cm)
Overall Length: .....	83" (210.8 cm)
Overall Width — Siderails Up: .....	36" (91.4 cm)
Overall Width — Siderails Stored: .....	30.25" (76.8 cm)
Siderail Length: .....	47" (119.4 cm)
Siderail Height Above Surface: .....	11" (27.9 cm)
Mattress Size: .....	26" x 75" (66.1 x 190.5 cm)
Maximum Head Elevation — Staff Controlled: .....	90°
Maximum Trendelenburg/Reverse Trendelenburg: .....	18°
Floor to Base Clearance: .....	3.5" (8.9 cm)
Caster Size: .....	8" (20.3 cm)
Maximum Weight Limit: .....	500 lbs. (226.8 kg)

**Direct Customer Service Hotline: USA 800-445-3730 | Canada 800-267-2337**

## Features

- Ideal for use in Trauma and Emergency
- X-ray capabilities from head-to-toe without moving the patient, saving time and reducing injury risk
- Full-length Lift Platform is spring loaded to accept all standard-size X-ray cassettes and grids. Lift Handles located on both sides can be activated to raise a cassette up to the sleep surface for optimal image quality
- Standard upright chest X-ray cassette holder and optional lateral X-ray cassette holder allow for increased efficiency
- Drainage Bags/Chest Tube Section mounts at all four corners
- TranStar® Stretchers feature the GentleRide™ shock-absorbing frame and siderail system for improved patient comfort and satisfaction.

## Accessories

- IV Transporter
- Patient Tray
- Restraint Straps
- Oxygen Tank Holder
- Liquid Oxygen Tank Holder
- Infusion Support System
- Siderail Covers
- Armboard
- Detachable Footboard
- Mattress Pad
- Lateral X-ray Cassette Holder
- Chart Holder

# Hill-Rom®

A HILLENBRAND INDUSTRY

Batesville, IN 47006 USA • 812-934-7777 • Fax 812-934-8189

Canada • 905-206-1355 • Fax 905-206-0561

International • 812-934-8173 • Fax 812-934-7191

www.hill-rom.com

CTG467



Patent Pending

Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2002 Hill-Rom Services, Inc.  
ALL RIGHTS RESERVED.