

TranStar® *Electric Stretcher*



The comfort and convenience of a hospital bed with the mobility of a stretcher

Siderail Controls allow patients to adjust stretcher for their individual comfort;

Staff Control Panel allows greater caregiver control

Steering Plus™ System for superior maneuverability and sharp cornering

OneStep® TuckAway Siderails for zero-gap transfers

Ergonomic Push Handles to reduce caregiver strain during transport

Hill-Rom®
A HILLENBRAND INDUSTRY

TranStar® Electric Stretcher

The Hill-Rom TranStar® Electric Stretcher is a highly mobile and durable unit that offers the patient the comfort and convenience of an electrical hospital bed while providing the mobility and features of the TranStar® line.

The TranStar® Electric Stretcher delivers new levels of mobility, comfort, and efficiency. It is the ideal solution for the extended-stay patient throughout the health care facility.

Electric Stretchers Value Package Configuration

Model	30" Width	Ergonomic Push Handles	IV Pole(s)	5th Wheel Steering	Knee Gatch	Mattress Options		
						4" Thick	5" Thick	30" Width
EL-550 Electric	x	x	1	x	x	Std	Option	Std

Labeling: **Standard decals include:** Emergency, Endoscopy, Cath Lab, GI Lab, PACU, Transport, Maternity, Radiology, Surgery, and Ambulatory Surgery

Custom decals: See your Hill-Rom Customer Service Agent or Account/Territory Manager if custom decals are desired.

Technical Data

Low Position — Floor to Top of Surface:	21.5" (54.6 cm)
High Position — Floor to Top of Surface:	35" (88.9 cm)
Overall Length:	83" (210.8 cm)
Overall Width — Siderails Up:	36" (91.4 cm)
Overall Width — Siderails Stored:	30.25" (76.8 cm)
Siderail Length:	47" (119.4 cm)
Siderail Height Above Surface:	14.5" (36.8 cm)
Mattress Size:	30" x 75" (76.2 x 190.5 cm)
Maximum Head Elevation — Staff Controlled:	90°
Maximum Head Elevation — Patient Controlled:	65°
Maximum Knee Elevation:	25°
Maximum Trendelenburg/Reverse Trendelenburg:	18°
Floor to Base Clearance:	3.5" (8.9 cm)
Caster Size:	8" (20.3 cm)
Volts:	120
Amperes:	4.75
Current Leakage:	<100 MicroAmps
Frequency:	50/60 Hz

Direct Customer Service Hotline: USA 800-445-3730 | Canada 800-267-2337

Hill-Rom®

A HILLENBRAND INDUSTRY

Batesville, IN 47006 USA • 812-934-7777 • Fax 812-934-8189

Canada • 905-206-1355 • Fax 905-206-0561

International • 812-934-8173 • Fax 812-934-7191

www.hill-rom.com

CTG468



Patent Pending

Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2002 Hill-Rom Services, Inc.
ALL RIGHTS RESERVED.

Features

- Ideal for use in PACU, ER Observation, and Outpatient
- Easy-to-read-and-activate Siderail Controls raise and lower the head and knee sections and allow for patient control for individual comfort
- Staff Control Panel, located at the footend, provides lockouts for the Siderail Controls and a control to raise the head section up to 90°. The CPR Levers, located on both sides of the unit, instantly lower the head section with a simple pull
- Fold-away Manual Cranks provide easy positioning when electrical power is not available
- TranStar® Stretchers feature the GentleRide™ shock-absorbing frame and siderail system for improved patient comfort and satisfaction.

Accessories

- IV Transporter
- Patient Tray
- Restraint Straps
- Oxygen Tank Holder
- Infusion Support System
- Siderail Covers
- Siderail Pads
- Detachable Footboard
- Convertible Footboard
- Mattress Pad
- Armboard
- Liquid Oxygen Tank Holder
- Lateral X-ray Cassette Holder
- Ankle Stirrups
- Chart Holder