

# TranStar®

## *Surgical Stretcher*



Articulating Head Section and Wrist Rest options to suit a variety of surgeon preferences

Steering Plus™ System for superior maneuverability and sharp cornering

OneStep® TuckAway Siderails for zero-gap transfers

Ergonomic Push Handles to reduce caregiver strain during transport

**Hill-Rom®**  
A HILLENBRAND INDUSTRY

# TranStar® Surgical Stretcher

The Hill-Rom TranStar® Surgical Stretcher, designed with the patient and staff in mind, combines mobility with functionality. The unit has features such as an articulating head section with multiple styles of head pads, a choice of two wrist rests, and a gas delivery/drape support device. Also, the unit has PACU Extenders that convert into armboards, allowing the caregiver to customize the stretcher to individual preferences.

These features, combined with a full array of accessories, make the TranStar® Surgical Stretcher the ideal solution for the operating room environment.

## Surgical Stretchers Value Package Configuration

Model	26" Width	Ergonomic Push Handles	IV Pole(s)	5th Wheel Steering	Knee Gatch	Mattress Options			
						4" Thick	26" Width	Flat Head Cushion	2 Piece Concave Head Cushion
SG-500 Surgical	x	x	1	x	x	Std	Std	Option	Option

Labeling: **Standard decals include:** Emergency, Endoscopy, Cath Lab, GI Lab, PACU, Transport, Maternity, Radiology, Surgery, and Ambulatory Surgery

**Custom decals:** See your Hill-Rom Customer Service Agent or Account/Territory Manager if custom decals are desired.

## Technical Data

Low Position — Floor to Top of Surface:	21.65" (54.9 cm)
High Position — Floor to Top of Surface:	34" (86.4 cm)
Overall Length:	92" (233.7 cm)
Overall Width — Siderails Up:	32" (81.3 cm)
Overall Width — Siderails Stored:	26.25" (66.7 cm)
Siderail Length:	47" (119.4 cm)
Siderail Height Above Surface:	14.5" (36.8 cm)
Mattress Size:	26" x 78" (66.0 x 198.1 cm)
Maximum Head Elevation:	70°
Maximum Knee Elevation:	25°
Maximum Weight Limit:	500 lbs. (226.8 kg)
Maximum Trendelenburg/Reverse Trendelenburg:	18°
Floor to Base Clearance:	3.5" (8.9 cm)
Caster Size:	6" (15.2 cm)

## Accessories

• Head Restraint Strap	• Mattress Pad	• Patient Tray
• Restraint Straps	• Oxygen Tank Holder	• Liquid Oxygen Tank Holder
• Foot Extender	• Chart Holder	• Convertible Footboard
• Infusion Support System	• Siderail Covers	• Armboard
• Lateral X-ray Cassette Holder	• Temporal Wrist Rest	• Integrated Removable Armboards and Pads
• Superior Wrist Rest	• Gas Delivery/Drape Support	• Hand Table

**Direct Customer Service Hotline: USA 800-445-3730 | Canada 800-267-2337**

# Hill-Rom®

A HILLENBRAND INDUSTRY

Batesville, IN 47006 USA • 812-934-7777 • Fax 812-934-8189

Canada • 905-206-1355 • Fax 905-206-0561

International • 812-934-8173 • Fax 812-934-7191

www.hill-rom.com

CTG469



Patent Pending

Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2002 Hill-Rom Services, Inc.  
ALL RIGHTS RESERVED.

## Features

- Ideal for use in Outpatient Surgery, OR, and Freestanding Surgery Centers
- Articulating Head Section and Wrist Rest options to suit a variety of surgeon preferences
- Unique Brake/Steer Pedal can "flip-up" off of the floor, providing maximum foot space for the physician as well as more room for equipment controls/foot pedals
- The Surgical Rail, mounted on both sides of the stretcher, provides greater flexibility by allowing for attachment of Gas Delivery/Drape Support Devices, PACU Extenders/ Armboards and Surgical Hand Tables
- TranStar® Stretchers feature the GentleRide™ shock-absorbing frame and siderail system for improved patient comfort and satisfaction.