

Hill-Rom Century™+ Bed



At a glance

- Durable solution meeting the flexibility of the medical/surgical environment.
- TuckAway siderails with embedded controls improve caregiver efficiency.
- Optional battery back-up provides ability to adjust bed during transport and power outages.
- Meets or exceeds the proposed FDA-sponsored Hospital Bed Safety Work Group (HBSW) dimensional guidelines for hospital beds.

Hill-Rom's Century™+ Bed is the ideal solution to the flexibility needs of today's health care providers. Manufactured to exact ISO 9001 quality standards and features TuckAway siderails with embedded controls. The bed delivers on the demands of the contemporary patient environment.

Starting from a feature-packed base model, caregivers can select options that significantly improve both staff efficiency and patient care quality. Caregivers can be assured the Century™+ Bed will fit into their budget constraints, and its distinctive features will reflect the quality of care they strive to provide for their patients.

Hill-Rom Century™+ Bed



Features:

- Nurse and patient siderail controls on both sides for ease of use.
- TuckAway siderails for easier patient transfers, easier linen changes and increased maneuverability.
- Siderail Hi-Lo control to make patient entry and exit safer and less strenuous.
- Trendelenburg and Reverse Trendelenburg controls with easy programming and engagement to allow the caregiver to be at the patient's side during special situations.
- Head elevation gauge to facilitate patient management.
- Patient lockout controls to deny patient access to any or all operating functions.

Options:

- Sidecom® UPG (upgradeable); NEL (Nurse Call Entertainment and Lighting); NB (Nurse and Bed Exit)
- Battery Backup
- Head & Footboards: Standard Century BLM (Blow Molded), available in Light Neutral Only.
- Laminated Options for Panels: Light Oak, Royal Oak, Bleached Oak, Natural Maple, Wild Cherry

Century™+ Bed Value Package Configurations

Model	Dust Cover	Auto Contour	5 UBS (5" urethane casters with central brake and steer)	5STL (5" santoprene total lock high mobility)	4SBS (4" santoprene with central brake and steer)	Night Light	CPR
CP 400	X	X	X			X	X
CP 300		X		X			X
CP 200					X		
CP 100				X			

Specifications:

Sleep Surface	.36" x 80" (91.4 x 203.5 cm)
Overall Frame (without bumpers)	.36" x 88" (91.4 x 223.5 cm)
Overall Length (with head/foot panels)	.88.5" (224.79 cm)
Overall Length (with head/foot panels and roller bumpers)	
Face Mount	.92.5" (234.95 cm)
Bottom Mount	.90.5" (229.87 cm)
Bedframe Clearance (with standard 5" casters)	
Rectangular Framework	.7.5" (19.05 cm)
Central Brake & Steer Cross Channel	.6.75" (17.14 cm)
Lift Arm Low Position	.2.875" (7.06 cm)
Overall Weight of Bed with blow molded head/foot panel, IV pole, and roller bumpers (without mattress)	.335 lbs. (151.95 kg)
High Position – Floor to Sleep Deck	.29.875" (75.64 cm)
Low Position – Floor to Sleep Deck	.18.125" (46.17 cm)
Maximum Head Elevation	.60°
Maximum Trendelenburg/Reverse Trendelenburg	.+/- 11°
Caster Size	.5" (12.7 cm)
Maximum Patient Weight Limit	.500 lbs. (226.8 kg)
Battery Back-Up (high-low cycles)	.100

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2003 Hill Rom Services, Inc.
ALL RIGHTS RESERVED.

CTG473ra 9/01/03