

Newborn Bassinet

PRODUCT SUMMARY



The bassinet is crafted from the highest quality materials and is designed for long lasting durability to create a home-like environment.

BASSINET
(P247)



Convenient handles on each side to facilitate transport of the neonate.



Two front locking casters to help keep bassinet securely in place.



Two shelves on either side of the bassinet are standard, providing easy access for documentation or placement of supplies such as formula bottles, diapers or linen.

NEWBORN BASSINET



SPECIFICATIONS:

Overall Cabinet and Handles with tub	39" W x 18" D x 40" H (99 cm x 46 cm x 102 cm)	Inside Drawer	26" W x 13" D x 4" H (66 cm x 33 cm x 10 cm)
Sleep Surface	26" W x 12" D x 2" H (66 cm x 30 cm x 5 cm)	Inside Door	27.5" W x 14" D x 13" H (70 cm x 36 cm x 33 cm)
Pull out Shelf	10" W x 13" D (includes handle) (25 cm x 33 cm)	Height of Mattress to Floor	33.5" / 85 cm
Chart Rack	11.5" W x 2" D x 8" H (29 cm x 5 cm x 20 cm)	Trend/Reverse Trend	10°

Hill-Rom
A HILLENBRAND INDUSTRY

Batesville, IN 47006 USA • 812-934-7777 • Fax 812-934-8189
Canada • 905-206-1355 • Fax 905-206-0561
International • 812-934-8173 • Fax 812-934-7191
www.hill-rom.com

CTG480

Due to limitations of the printing process, colors may vary slightly. When specifying, obtain color samples by contacting the fabric or finish manufacturer directly.

Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
©2002 Hill-Rom Services, Inc.
ALL RIGHTS RESERVED.