

Hill-Rom "Flex" Highback and Highback Chairs



The "Flex" Highback and Highback Chairs are designed for sitting comfort.



Straight Leg Highback "Flex" Chair



Bentwood Highback "Flex" Chair

At a glance

- All chairs feature the ease-of-cleaning and quality of construction that have made Hill-Rom chairs the long-standing favorite of caregivers throughout the world.
- "Wall saver" design.
- Designed for comfort and support.
- Choose traditional Straight Leg or Bentwood for added style and flare.
- Compactly designed for today's patient rooms.



Bentwood Highback Chair



Straight Leg Highback Chair

Hill-Rom "Flex" Highback and Highback Chairs



Specifications:

	Bentwood	Straight Leg	Bentwood "Flex"	Straight Leg "Flex"
Floor Print	24" W x 30" D (61 cm x 76 cm)	23" W x 30" D (59 cm x 76 cm)	24" W x 30" D (61 cm x 76 cm)	23" W x 30" D (59 cm x 76 cm)
Overall Height	41.5" (106 cm)	41.5" (106 cm)	41.5" (106 cm)	41.5" (106 cm)
Overall Width	24" (61 cm)	23" (59 cm)	24" (61 cm)	23" (59 cm)
Overall Depth	30" (76 cm)	30" (76 cm)	30" (76 cm)	30" (76 cm)
Seat Height	18" (46 cm)	18" (46 cm)	18" (46 cm)	18" (46 cm)
Seat Width	19" (48 cm)	20" (51 cm)	19" (48 cm)	20" (51 cm)
Seat Depth	18" (46 cm)	18" (46 cm)	18" (46 cm)	18" (46 cm)
Arm Height	25" (64 cm)	24" (61 cm)	25" (64 cm)	24" (61 cm)
Shipping Weight	45 lbs (21 kg)	45 lbs (21 kg)	55 lbs (25 kg)	55 lbs (25 kg)
COM Yards	4 yds	4 yds	4 yds	4 yds

Stains




Natural Maple 31M



Light Oak R30



Medium Oak 820



Dark Cherry H4C



Wild Cherry T9C

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2003 Hill Rom Services, Inc.
ALL RIGHTS RESERVED.

CTG481ra 10/10/03