

Recliner Chairs - Bentwood and Straight Leg

PRODUCT SUMMARY



The Recliner Chair is designed for comfort, stability and efficiency.

BENTWOOD RECLINER

(REC-550-31m-604)

STRAIGHT LEG RECLINER (SHOWN BELOW)

(REC-500-31m-604)



3 specific chair positions



Easy grip armrests stabilize chair during entry and exit.



Front locking casters keep chair stationary during entry and egress.




Easy to replace cushions for convenient cleaning and servicing.



Open seat/back design for comfortable seating and cleaning.



Food tray available as an accessory



SPECIFICATIONS:

Overall Height	46" / 117 cm
Overall Width	29" / 74 cm
Overall Depth	33" / 84 cm
Shipping Weight	94 lbs / 43 kg
Seat Height	19" / 48 cm
Seat Width	24" / 61 cm
Seat Depth	18" / 46 cm
Arm Height	25" / 64 cm
Floor Print	29" W x 33" D (74 cm x 84 cm)
Extended	29" W x 59" D (74 cm x 150 cm)

Hill-Rom
A HILLENBRAND INDUSTRY

Batesville, IN 47006 USA • 812-934-7777 • Fax 812-934-8189
Canada • 905-206-1355 • Fax 905-206-0561
International • 812-934-8173 • Fax 812-934-7191
www.hill-rom.com

Due to limitations of the printing process, colors may vary slightly. When specifying, obtain color samples by contacting the fabric or finish manufacturer directly.

Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
©2002 Hill-Rom Services, Inc.
ALL RIGHTS RESERVED.