

Procedure Stool

PRODUCT SUMMARY



Designed to adjust from a low of 20" to 24" for convenience and comfort, using a foot-controlled gas spring. A five caster base, with two locking casters, provides fluid mobility.

PROCEDURE STOOL
(P118B524)

The Procedure Stool has the following standard features:



Two locking casters



Soft, wide vinyl seat



Foot controlled gas spring adjustment

SPECIFICATIONS:

Seat Diameter	15" / 38 cm
Seat Travel Distance	4" / 10 cm

Hill-Rom[®]
A HILLENBRAND INDUSTRY

Batesville, IN 47006 USA • 812-934-7777 • Fax 812-934-8189
Canada • 905-206-1355 • Fax 905-206-0561
International • 812-934-8173 • Fax 812-934-7191
www.hill-rom.com

Due to limitations of the printing process, colors may vary slightly. When specifying, obtain color samples by contacting the fabric or finish manufacturer directly.

Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2002 Hill-Rom Services, Inc.
ALL RIGHTS RESERVED.