

TranStar[®]

Value Procedural Stretcher

(Models PC-100/PC-200)



Easy access four-corner brake and steer pedals
Ergonomic push handles* help reduce caregiver strain during transport
Highly mobile, economical and durable patient transport

** Standard on PC-200 model procedural stretchers*

Hill-Rom[®]
A HILLENBRAND INDUSTRY

TranStar[®]

Value Procedural Stretcher

(Models PC-100/PC-200)

The PC-100 and PC-200 TranStar[®] Procedural Stretchers are Hill-Rom's most economical stretchers ever. Maneuverability is improved as this lightweight stretcher offers four-corner brake/steer controls that are accessible from any position around the stretcher. The ergonomic push handles (PC-200 package) help reduce caregiver strain during transport. The siderails are designed to minimize the opportunity for pinched fingers and tubes. In critical situations, the stretcher can also be easily placed into emergency Trendelenburg position by simply lifting the foot end of the stretcher.

When the goal is value and economics in a maneuverable, efficient and durable stretcher, the TranStar PC-100 and PC-200 packages are the clear answer.

Procedural Stretchers PC-100 and PC-200 Value Package Configurations

Model	26" Width	4-corner Brake & Steer Controls	Ergonomic Push Handles	IV Pole(s)	Mattress Options				
					3" Thick	4" Thick	5" Thick	26" Width	30" Width
PC-200	x	x	x	1	Std	Op	Op	Std	
PC-100	x	x			Std	Op	Op	Std	

Note: CAL 129 Flammability requirement is an option on all mattresses.

Labeling: **Standard decals include:** Emergency, Endoscopy, Cath Lab, GI Lab, PACU, Transport, Maternity, Radiology, Surgery, and Ambulatory Surgery.

Custom decals: See your Hill-Rom Customer Service Agent or Account/Territory Manager if custom decals are desired.

Technical Data PC-100—PC-200

Low Position — Floor to Top of Surface:	23" (58.4 cm)
High Position — Floor to Top of Surface:	35" (88.9 cm)
Overall Length:	83" (210.8 cm)
Overall Width:	32" (81.3 cm)
Siderail Length:	47" (119.4 cm)
Siderail Height Above Mattress/Surface:	14" (36.8 cm)
Mattress Size:	26" (66 cm) x 75" (190.5 cm)
Maximum Head Elevation:	90°
Maximum Trendelenburg/Reverse Trendelenburg:	18°
Floor to Base Clearance:	3.5" (8.9 cm)
Caster Size:	8" (20.3 cm)
Maximum Weight Limit:	500 lbs. (226.8 kg)

Direct Customer Service Hotline: USA 800-445-3730 | Canada 800-267-2337

Hill-Rom[®]

A HILLENBRAND INDUSTRY

Batesville, IN 47006 USA • 812-934-7777 • Fax 812-934-8189

Canada • 905-206-1355 • Fax 905-206-0561

International • 812-934-8173 • Fax 812-934-7191

www.hill-rom.com

CTG537



Patent Pending

Hill-Rom reserves the right to make changes without notice in design, specifications, and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2002 Hill-Rom Services, Inc.
ALL RIGHTS RESERVED.

Features

- The economic choice for Emergency, OR, PACU, and Outpatient
- Four-corner brake and steer controls allow easy access from any position
- Low height improves on/off transfers and increases patient safety
- Ergonomically correct transport handles help reduce caregiver strain during patient transport (PC-200 model)
- TranStar[®] Stretchers features the GentleRide™ shock-absorbing frame and siderail system for improved patient comfort and satisfaction

Accessories

- IV Transporter
- Patient Tray
- Utility Tray
- Restraint Straps
- Oxygen Tank Holder
- Liquid Oxygen Tank Holder
- Infusion Support System
- Siderail Covers
- Armboard
- Detachable Footboard
- Mattress Pad
- Lateral X-ray Cassette Holder
- Chart Holder
- Convertible Footboard