

**Hill-Rom®**

Art of Care® Bedside Cabinets



**Metropolitan**

Product **BCSM** Base Finish **31**  
Laminate **734** Powder Coat **QSL** Casters **C**

## Hill-Rom® Bedside Cabinets: A marriage of beauty and durability.

Built for the specific needs of healthcare users, Hill-Rom's Metropolitan™, Artisan™, and Aero™ bedside cabinets deliver a striking combination of function and form. Hill-Rom's Metropolitan and Aero bedside cabinets are constructed around a durable steel frame, while Hill-Rom's Artisan cabinet is made with an all-wood construction and steel sub-bracketry for structural strength. Each cabinet has a high pressure laminate top and plastic inner drawers for easy cleaning. Each Hill-Rom bedside cabinet has the revolutionary Healthcoat 222™ process that bakes and seals the finish directly to the wood, making it able to withstand years of repeated cleaning. Even the hardware promotes easy cleaning without compromising the aesthetic appeal. The inherent function, beauty, and reliability of Hill-Rom's Metropolitan, Artisan, and Aero bedside cabinets add lasting warmth and style to your environment.



**Aero**

Product **BSCR** Base Finish **HM**  
Laminate **WM9** Casters **C**



**Artisan**

Product **BSCA** Base Finish **CC**  
Laminate **117** Casters **C**

### At a glance

- Three unique aesthetic styles
- Hardwood and hardwood veneer drawer fronts and side panels
- High pressure laminate tops for easy cleaning
- Hill-Rom® Healthcoat 222™ Finish for Healthcare environments

**Hill-Rom®**

Enhancing Outcomes for Patients and Their Caregivers.™

# Hill-Rom®

## Art of Care® Bedside Cabinets



### Easy Clean Hardware

Carefully selected to provide beauty and performance.



### Generously Rounded Corners in Storage Drawers

Eliminates corner crevices for easy cleaning.



### Easy-Roller Casters

Caster design makes moving cabinets simple.



### Screw-In Glides


Our optional world-class glides stand up to heavy use and allow for adjustment.

## Technical Specifications

### BSCM

Metropolitan™ Cabinet

W 25 D 19 H 32.5



### BSCA

Artisan™ Cabinet

W 20.5 D 18.5 H 32.5




### BSCR

Aero™ Cabinet

W 21.5 D 17.5 H 30.5

Also matches Opera Style



### Product

### Wood Finish

### Laminate

### Powder Coat

### Options

Please Specify the Following for Each BEDSIDE CABINET

Please specify Product from list above by Four Digit Product Code.  
For any product not shown above, please submit a Request for Special.

See Hill-Rom® Finish Options brochure<sup>1</sup> or visit [www.hill-romdesign.com](http://www.hill-romdesign.com)

See Hill-Rom® Finish Options brochure<sup>1</sup> or visit [www.hill-romdesign.com](http://www.hill-romdesign.com)

Metropolitan only – See Hill-Rom® Finish Options brochure<sup>1</sup> or visit [www.hill-romdesign.com](http://www.hill-romdesign.com)

Casters Standard. Glides available for Metropolitan and Artisan only.

<sup>1</sup>CTG255rb rev 3

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2008 Hill-Rom Services, Inc ALL RIGHTS RESERVED.

CTG731 rev 2 5/16/08



Enhancing Outcomes for Patients and Their Caregivers.™

USA 800-638-2546  
Canada 800-267-2337

[www.hill-rom.com](http://www.hill-rom.com)