

Hill-Rom Contour™ Headwall System



A headwall unit that meets your basic service needs in a trim, effective way

Have an environment that doesn't require the functionality and flexibility that Hill-Rom's other great headwalls offer? Contour meets your minimum functional needs for the right price. You will also be amazed by its sleek design, ease of installation and high quality.

Aesthetics

- Recessed to take focus off of the headwall and place it where it belongs, on the patient
- Less impact on the space of the room
- Choose from a variety of laminate colors to meet your room's aesthetic needs

Price

- Overall price competitive with inwall

At a glance

- Offers a sleek, trim look that minimizes the hospital feel of the room
- Priced to fit your budget
- Recessed into wall to limit impact on room space
- Tested to ensure it meets safety and quality standards
- Can install before or after drywall, thus offering greater flexibility and reduced opportunity for damage to product
- Behind schedule? We have existing stock of already-made, preconfigured units that can be shipped within one to two weeks of receiving the order

Hill-Rom Contour™ Headwall System



Function

- Equipment management rails allow you to manage your basic equipment needs
- Installation is quick and easier than traditional inwall
- Need more services? Contour™ headwall system comes in a longer version that offers four additional duplexes and one additional low voltage provision

Quality

Safety and quality testing ensures product meets or exceeds current standards without additional cost.

▶ Specifications

	Standard	Stock	Extended
Gases (oxygen, air or vacuum)	Up to 6 (configurable)	3 (1 Oxygen, 1 Air, 1 Vacuum)	Up to 6 (configurable)
Electric (in # of duplexes)	2 normal, 2 emergency	2 normal, 2 emergency	4 normal, 4 emergency
Low Voltage Provisions	3	3	4
Nurse Call	1, 2, 3 or 4 gang	4 gang	1, 2, 3 or 4 gang
Length	53.5"	53.5"	59.5"
Width	16.625"	16.625"	16.625"
Distance from Floor to Top of Unit	79.25"	79.25"	79.25"
Distance from Floor to Bottom of Unit	25.75"	25.75"	19.75"
Distance from Floor to Center Line of Top Row of Gases	60"	60"	60"
Wall Depth	3.5" deep studs	3.5" deep studs	3.5" deep studs
Wall Depth (Back to Back)	6" deep studs	6" deep studs	6" deep studs
Lead Time	Normal	Expedited (1 – 2 weeks)	Normal

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
©2004 Hill-Rom Services, Inc.
ALL RIGHTS RESERVED.
CTG733ra 7/27/04

